



ANALIZA BIAŁEK



100008784

## Laboratorium Analizy

**Client:** Emperor  
**Order number:** 100008784  
**Name:** EMPEROR GW-501516  
**Number of samples:** 1 container

**Batch number:** -----  
**Receive date:** 2025-07-11  
**Start date:** 2025-07-25  
**End date:** 2025-07-28

## Sample identification

Substance	Method of taking a sample	Form
GW-501516	Sample provided by customer	White capsules in a plastic container

## Results

Substance	Sample name	Content	Purity
GW-501516	EMPEROR GW-501516 10MG/CAPS	11,02mg/caps	-----

## Research comment

The test sample was applied in the amount of 2 $\mu$ l to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H<sub>2</sub>O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peak of the sought substance appeared. Its concentration, when compared to the reference chromatogram, is 11,02mg/tab. Height the measurement error is + - 0,05 mg. The actual state of the sample sent before the test is presented on the photo below

## Photo of the test sample



## Password

2N28S7G3

**Analiza Białek sp. z o.o.**  
**KRS:** 0000966168  
**NIP:** 8992922970  
**REGON:** 521751875

**Headquarters and Laboratory:**  
ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**  
Analysis and research specialist  
mgr Arkadiusz Zajac




ANALIZA BIAŁEK



100008784

## Chromatogram



## Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	95%	5%
2min	95%	5%
12min	5%	95%
15min	5%	95%

## Password

2N28S7G3

**Analiza Białek sp. z o.o.**

**KRS:** 0000966168

**NIP:** 8992922970

**REGON:** 521751875

**Headquarters and Laboratory:**

ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**

Analysis and research specialist  
mgr Arkadiusz Zajac